



MK-dent



CELEBRATING
25 YEARS

PROMOTION

NOV 20 - JAN 21



FAIRCARE

WARRANTY

- FairCare Warranty is an **Unconditional guarantee** which includes all parts and labour during the warranty period (except physical damage) which is unique within the Australian and New Zealand market.
- FairCare Warranty provides a swift and streamlined service to get your instrument back in action, through the use of genuine MK-dent replacement parts.
- FairCare offers you competitive prices on genuine replacement parts for post warranty repairs, ensuring the safety and longevity of your instrument.

Take further advantage of the discounted prices inside with this **AMAZING** offer



P R O M O T I O N 2 0 2 0

INTEREST FREE¹

When you spend from this catalogue

\$ 5,500 to \$ 6,500 - 6 months 0 %

\$ 6,501 to \$ 7,500 - 9 months 0 %

\$ 7,501 onwards - 12 months 0 %

with equal monthly payments

- ✓ Instant approval in most cases
- ✓ No finance setup cost
- ✓ Often no forms to complete
- ✓ Easy cash-flow solution

¹In relation to any promotions or offers by, relating to or referencing Credabl Pty Ltd (ACN 615 968 100) (Credabl) that are contained in this document (Credabl Promotion), this document is a guide only and does not constitute any recommendation on behalf of Credabl. The information in this document is general in nature and Credabl has not taken into account your personal objectives or financial circumstances or needs when preparing it. Before acting on this information you should consider if it is suitable for your personal circumstances. Credabl is not offering financial, tax or legal advice. You should obtain independent financial, tax or legal advice as appropriate. Credabl reserves the right to cease offering these products at any time without notice. The Credabl Promotion is subject to Credabl's terms and conditions, fees and charges and eligibility criteria. MK-Dent Australia Pty Ltd ACN 161 945 946 (MK-Dent) is not an agent or partner of Credabl and is a separate and independent entity from Credabl. In no event will Credabl be liable for any actions, representations and / or decisions of MK-dent or its related bodies corporate representatives, directors, officers, employees, agents or advisors. The issuer and credit provider of the Credabl Promotion is Credabl Pty Ltd (ACN 615 968 100) Australian Credit Licence No. 499547.



PRIME LINE



AeroAcoustic blade



Refined ceramic bearings



Automatic pressure regulation¹



Reverse flow protection



Glass rod fiber optics



Scratch resistance
Iconic Prime body

25
W

Powerful & Precise
Power Head - 25W

55
dB

Incredibly Quiet
Only 55 dB noise



0%
FINANCE
AVAILABLE

Details
on page 2



Now with 2 Years
FairCare Warranty²
Unconditional warranty
except physical damage



Coupling Connection	Power Head 25 Watt	Small Head 18 Watt
Kavo	HP21KL	HP22KL
NSK	HP21NL	HP22NL
W&H	HP21WL	HP22WL
Sirona	HP21SL	- N/A -

RRP: \$ 1,595.00

- Buy 1 x Prime Line - Get 25% OFF
Pay only \$1,196.25 - SAVE \$398.75
- Buy 3 x Prime Lines - Get 30% OFF
Pay only \$1,089.00ea - SAVE \$1,520.00

1) Only available for Kavo coupling

2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

ECO LINE TURBINES



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Push Button system
- Smooth Eco body

- \$847.00 (Fibre Optic)
- \$627.00 (Non Optic)



BUY 3
\$2,145
Save \$396

Fibre Optic



Fibre Optic

LB11/11

BUY 4
GET 2 FREE
Slow Speeds
(LB11/11)
\$3,245
Save \$913



Fibre Optic

LB11/11

BUY 5
GET 3 FREE
Slow Speeds
(LB11/11)
\$4,235
Save \$1,155

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)
Kavo	HE21KL	HE22KL
NSK	HE21NL	HE22NL
W&H	HE21WL	HE22WL
Sirona	HE21SL	- N/A -

ECO LINE CONTRA - ANGLES



1 Year FairCare Warranty²
Unconditional warranty
except physical damage on
ECO LINE turbine & contra-angle

\$ Buy 3 - Get 15% Off Mix & Match



SAFE-TECH PUSH BUTTON



- 1:5 Red ● Max. 200,000 rpm Internal water
- 1:1 Blue ● Max. 40,000 rpm Internal water
- 5,4:1 Green ● Max. 7,407 rpm Internal water
- 1:1 Blue ● Max. 40,000 rpm Internal water

- LE15L (F/Optic) RRP: \$1,485.00 **Buy 3 - \$1,262.25 ea**
- LE15 (Non-optic) RRP: \$1,287.00 **Buy 3 - \$1,093.95 ea**
- LE11L (F/Optic) RRP: \$984.50 **Buy 3 - \$836.80 ea**
- LE11 (Non-optic) RRP: \$759.00 **Buy 3 - \$645.15 ea**
- LE54LP (F/Optic) RRP: \$1,155.00 **Buy 3 - \$981.75 ea**
- LE54P (Non-optic) RRP: \$880.00 **Buy 3 - \$748.00 ea**
- LE01L (F/Optic) RRP: \$935.00 **Buy 3 - \$794.75 ea**
- LE01 (Non-optic) RRP: \$748.00 **Buy 3 - \$635.80 ea**

1) Eco Line Power Head 24 Watt 2) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019



**Buy More - Save More
Mix & Match**



**1 year Standard warranty on Basic Line
Slow-speeds, Airmotors & Couplings**



**NEW
MODEL**




-  1:1 Blue ●
-  Max. 20,000 rpm
-  Push button

● **LB11/11**

Buy 1 - \$385.00

Buy 3 - \$329.00 ea






-  4:1 Green ●
-  Max. 5,000 rpm
-  Snap-on Cups

● **LB41/11P**

Buy 1 - \$462.00

Buy 3 - \$396.00 ea



-  4:1 Green ●
-  Max. 5,000 rpm
-  Push button

● **LB41/11**

Buy 1 - \$539.00

Buy 3 - \$495.00 ea



-  1:1 Blue ●
-  Max. 20,000 rpm

● **LB02**



Buy 1 - \$385.00

Buy 3 - \$329.00 ea

AIRMOTORS 



**Shorter
and
Lighter**

-  Internal water
-  LED

● **AM1006**
\$ 990.00



-  Internal water

● **AM1004**
\$ 715.00



-  External water

● **AM0014**
\$ 429.00



● **AM0004**
\$ 429.00

COUPLINGS 



-  MK-dent & Kavo
with LED & water

● **QC6016KW**
\$ 495.00



-  NSK connection
with LED

● **QC6016NT**
\$ 352.00



-  Sirona connection
with LED & water

● **QC6016SW**
\$ 495.00



-  W&H connection
with optic

● **QC5016W**
\$ 385.00

▼ PROPHY LINE ULTRASONIC SCALER



**Buy 3 - Get 15% OFF
Plus
FREE tip with every scaler**



1 year warranty



EMS Non-optic type EN-041

- SC21E (SetSC2)
Buy 1 - \$429.00
Buy 3 - \$364.65 ea



EMS Piezon Non-optic type EN-061

- SC21E2 (SetSC5)
Buy 1 - \$528.00
Buy 3 - \$448.80 ea



EMS Piezon LED type EN-060A

- SC21EL (SetSC1)
Buy 1 - \$682.00
Buy 3 - \$579.70 ea

NEW



KaVo Piezo LED type 1.007.3995

- SC21KL* (SetSC6)
Buy 1 - \$682.00
Buy 3 - \$579.70 ea



Acteon-Satalec Non-optic type F 12 281

- SC21S (SetSC4)
Buy 1 - \$429.00
Buy 3 - \$364.65 ea



Acteon-Satalec LED type F 12 609

- SC21SL (SetSC3)
Buy 1 - \$682.00
Buy 3 - \$579.70 ea

*Only for use with MK-dent STxxE/ED or EMS original scaler tips

▼ PROPHY LINE AIR POLISHING DEVICE



Get 3 bottles of Prophy Powder FREE with every Unit



1 year warranty

- PR1011
\$1,287.00



Prophy Powder (300 gms)
\$ 60.50 each

- Mint - PP1011
- Lemon - PP1012
- Strawberry - PP1013

- Ergonomically shaped, compact and light-weight air-driven polishing device

- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation

Coupling Connection:

Kavo - PR1011K, NSK - PR1011N
W&H - PR1011W, Sirona - PR1011S



**Buy 4 - Get 10% OFF
Mix & Match**

- ✓ High-quality DLC coating
- ✓ Extremely hard and wear-resistant surface up to 40% longer service life



**SC21KEY
Torque Key**
\$ 49.50

• SCALING DLC



**ST01ED/ST01SD
SATELEC 1**

\$88.00



**ST02ED/ST02SD
EMS P
SATELEC 1S**

\$88.00



**ST03ED/ST03SD
EMS A**

\$88.00



**ST04ED/ST04SD
EMS C
SATELEC 2**

\$132.00

• PERIO DLC



**ST06ED/ST06SD
EMS PS
SATELEC 10Z**

\$132.00



**ST10ED/ST10SD
SATELEC 10P**

\$154.00



**ST08ED/ST08SD
SATELEC S07**

\$154.00

• ENDO DLC

• SCALING STAINLESS STEEL



**ST01E/ST01S
SATELEC 1**

\$63.80



**ST02E/ST02S
EMS P
SATELEC 1S**

\$63.80



**ST03E/ST03S
EMS A**

\$63.80



**ST04E/ST04S
EMS C
SATELEC 2**

\$88.00

• PERIO STAINLESS STEEL



**ST06E/ST06S
EMS PS
SATELEC 10Z**

\$99.00



**ST07E/ST07S
EMS DPL3
SATELEC 10X**

\$110.00



**ST09E/ST09S
EMS HPL**

\$110.00



**ST10E/ST10S
SATELEC 10P**

\$110.00



**ST11E/ST11S
EMS RT1**

\$110.00

• ENDO STAINLESS STEEL



**ST05E/ST05S
EMS D
SATELEC 5AE**

\$88.00



**ST08E/ST08S
SATELEC S07**

\$110.00



**ST12E/ST12S
120°
FILE HOLDER**

\$99.00



**ST13E/ST13S
90°
FILE HOLDER**

\$99.00



**ST14E/ST14S
180°
FILE HOLDER**

\$99.00



**ST15E/ST15S
SATELEC ET40**

\$110.00



**ST16E/ST16S
SATELEC ET40D**

\$110.00



**ST17E/ST17S
SATELEC ET20**

\$110.00



**ST18E/ST18S
SATELEC ET20D**

\$110.00

• CAVITY PREPARATION STAINLESS STEEL



**ST19E/ST19S
SATELEC EX1**

\$132.00



**ST20E/ST20S
EMS SBM
SATELEC EX2**

\$132.00



**ST21E/ST21S
EMS SBM
SATELEC EX3**

\$132.00



**ST22E/ST22S
SATELEC EXL**

\$132.00



**ST23E/ST23S
SATELEC EXR**

\$132.00

✔ Fully Synthetic Premium Service Oil



- ✔ For any brand of dental handpieces and motors
- LU1011 (500 ml)
\$ 59.40



- ✔ For Kavo Quattrocare 2 lubrication system
- LU1022 (500 ml)
\$ 59.40

MANUFACTURER

MK-dent

MK-dent GmbH
Marie-Curie-Str. 2
22941 Bargteheide
Germany

www.mk-dent.com
info@mk-dent.com

YOUR SPECIALIST DEALER



3/ 1-7 Union St. Stepney, SA 5069
Toll Free 1800 882 412
www.citydentalsupplies.com.au



CELEBRATING
25 YEARS



■ MK-DENT PARTNERS



Global network of authorized dealers



UNIQUE
FAIRCARE
WARRANTY



- MK-dent was founded in 1995 in Bargteheide, Germany and has developed into one of the premium dental instrument manufacturers.
- MK-dent products are currently available in more than 100 countries through an extensive distribution network.

• All catalogue prices are inclusive of GST Promotion valid until 31st Jan 2021 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.